

Inhibition Panel – Haemolytic Specimens

Intended Use

The haemolytic specimens are designed to perform standard interference testing by diagnostic manufacturers and researchers during assay development, evaluation and post-marked surveillance. The panel suited for execution of routine investigations for the survey of inhibition effects by haemoglobin in serologic/immunologic/biochemical assays as well as nucleic acid amplification techniques (NAT). This panel is for Research Use Only.

Product Description

The most frequently performed interference studies are for the following serum indices: haemolysis, lipaemia and icterus.¹ Beside haemolysis, lipaemia and increased bilirubin concentrations may also interfere with determinations, whereby bilirubin is the second most common interference after haemolysis.^{2,3} The haemolytic specimens consist of human plasma, which contain intracellular components of erythrocytes and therefore increased values of haemoglobin.

The haemolytic specimens are present in low, medium and high concentrations and a negative control. The four components originated from one normal healthy donor. The haemolytic plasma is a completely native product. It originates from K₂EDTA (1.8 mg/mL) plasma and contains a specific amount of haemoglobin of lysed erythrocytes from the same individual. The measurements were performed by use of Sysmex KX-21N. Our product mimics haemolysis in a very realistic way because the intracellular components are also included.⁴ Furthermore, no preservatives were added.

Article Code	Components	Concentration Range	Quantity
IHH-005	Haemolytic negative control	0.0 - 0.1 g/dL	each 1.0 mL
IHH-010	Haemolytic plasma low	0.5 - 0.7 g/dL	
or	Haemolytic plasma medium	0.9 - 1.1 g/dL	
IHH-XXX	Haemolytic plasma high	1.9 - 2.2 g/dL	

Panels with 5, 10 or more individuals (Article Code: IHH-005, IHH-010 or IHH-XXX) contain the corresponding quantity 5 or 10 x 1.0 mL resp. custom-made. All components are ready to use, if requested only the analyte has to be added into the specific matrix in requested concentration.

The Inhibition Panel provides great flexibility in usage in that any sample preparation platform and subsequent analytical assay may be validated. The volumes of the final samples spiked with analyte and/or the final concentration of analyte may be adjust by end user so that they are within requested levels for the assay of interest.

Storage and Stability

The Inhibition Panel is stable until the expiration date, provided it is unopened and stored at a minimum of -20°C. It is recommended to store the product at -20°C to -80°C to ensure highest quality. The specimens can be aliquoted in smaller volumes to avoid repeated freeze/thaw cycles. After usage, promptly refreeze to avoid leaving the product thawed for a longer period of time. A minimal loss of activity is possible. Ensure proper storage and limit exposure time.

Instructions for Use

The Inhibition Panel was designed for usage with any assay system, which detects analytes derived or stored in whole blood or blood derivatives. All specimens are ready to use and only analyte has to be added into the specific matrix in requested concentration. It is recommended that the analyte should always be added with an equal volume/concentration, so that

inhibition/non-inhibition effects (interference/non-interference claims) of the extraction and/or assay system are consistently verified.

The haemolytic specimens should be thawed at room temperature and extensively mixed by inverting to ensure homogenous suspensions. Avoid foaming and contamination of components when opening and aliquoting. After usage, return immediately to storage conditions to prevent long periods of thawing. The analyte has to be in a physiologic buffer, plasma, serum or other blood derivative to be compatible with the panel specimens.

The analyte can be added directly to the specimen vials, or the specimens can be aliquoted in smaller volumes and then spiked to a volume within the range of 1% to 50% of the total volume, depending on the concentration that is used for spiking. It should be noted that the initial concentration of components decreases accordingly. The added volume of the analyte should not exceed more than 50% of the total volume, as this will reduce the effectivity range of the potential interfering substance.

After the analyte is spiked, proceed with the samples as described by extraction or assay instructions from manufacturer. Follow the appropriate instructions and recommendations for safe handling and testing of human materials.

Limitations

The Inhibition Panel – haemolytic specimens are not intended to replace internal diluents, controls, standards or calibrators of *in vitro* diagnostic (IVD) test kits nor for blood bank screening. This product is for Research Use Only and not for human or animal diagnostics, or for therapeutic use.

Do not use this product beyond the expiration date. In the event of microbial contamination or high turbidity, do not use. If this product is used in tests, please follow the exact instructions of the manufacturer. Any deviations from the manufacturer instruction for use produce unreliable results.

Each laboratory has the responsibility to ascertain the suitability of the Inhibition Panel components for its particular application and to establish their own guidelines for interpretation of results. Data are provided for informational purpose. It is not claim that others can duplicate test results exactly.

Precautions and Disposal Note

Although whole source material has tested negative/non-reactive for anti-HIV, anti-HCV and HBsAg by CE-marked tests, they should be considered as potentially infectious, since no testing method can absolutely exclude the potential danger of infection. Treat this material with the same precaution as potentially infectious material.^{5,6}

During the processing of the samples, longer refrigeration times of one to three days may lead to an increase in, or reduction of certain substances by degradation processes. All material should be handled and discard in a safe manner that complies with local and national regulations and will inactivate pathogenic agents.

References

1. Dimeski G., Interference Testing. Clin Biochem Rev. 2008 Aug; 29 Suppl 1:S43-48.
2. Halwachs-Baumann G., Labormedizin Klinik-Praxis-Fallbeispiele, Springer-Verlag/Wien 2011, 2. Auflage
3. Nikolac N., Lipemia: causes, interference mechanisms, detection and management. Biochem Med 2014;24(1):57-67. Doi: 10.11613/BM.2014.008. eCollection 2014. Review.
4. Thomas L., Haemolysis as influence and interference factor. eJIFCC vol 13 no 4: <http://www.ifcc.org/ifccfiles/docs/130401002end.pdf>

5. Ausschuss für Biologische Arbeitsstoffe, Technische Regeln für Biologische Arbeitsstoffe TRBA 100: Schutzmaßnahmen für Tätigkeiten mit biologischen Arbeitsstoffen in Laboratorien; Ausgabe Oktober 2013, GMBI 2014, Nr. 51/52 vom 17.10.2013, Änderung vom 30.06.2014, GMBI Nr. 38, <http://www.baua.de/cae/servlet/contentblob/673098/publicationFile/48545/TRBA-100.pdf>
6. Richtlinie 2000/54/EG des Europäischen Parlaments und des Rates über den Schutz der Arbeitnehmer gegen Gefährdung durch biologische Arbeitsstoffe bei der Arbeit (Siebte Einzelrichtlinie im Sinne von Artikel 16 Absatz 1 der Richtlinie 89/391/EWG), vom 18. September 2000 (ABl. EU Nr. L 262 S. 21) http://www.gaa.baden-wuerttemberg.de/servlet/is/16050/2_1_07.pdf